## NCIMS Appendix N Pilot Program

For all NCIMS testing, refer to the Tetracycline Pilot Program Bulk Milk Tanker Screening Test form for IDEXX SNAP Tetracycline (Appendix N Pilot Program).

- Upon initial screen using the Procedure, a positive result indicates an Initial Positive.
- For milk samples yielding an Initial Positive result, promptly retest the SAME sample in duplicate using the sample dilution protocol below. Test a positive and a negative control along with the duplicate retest samples. Do not use checksets.
- To obtain detection levels reported in the "Performance Information" section of this insert, the sample must be diluted with previously screened tetracycline negative milk. In a clean container, mix 1 part of sample milk to be tested with 9 parts of negative milk. For example, if using calibrated pipettes, add 0.5 mL of sample milk to 4.5 mL of negative milk. Mix well and proceed to follow the test kit instructions under the "Procedure" section of this insert.
- The preparation of the positive and negative control samples is described in the "Checksets, Positive and Negative Controls in the United States for NCIMS Testing" section.

### Interpretation of Retest Results

If the controls test correctly, negative results for both retest samples indicate a Negative Result (report as Not Found).

If the controls test correctly and one or both of the retest samples test Positive, the sample is a "Tetracycline Presumptive Positive Test."

In accordance with Appendix N Pilot Program, report all presumptive positives to the appropriate regulatory agency.

If either control does not test correctly, then repeat the assay. If either control still does not test correctly, do not continue testing. Contact the appropriate regulatory agency to determine milk status. Contact IDEXX Technical Services for technical assistance.

### Performance Information

#### Sensitivity (using dilution protocol)

Dose Response Information	% Positive		
	Tetracycline	Oxytetracycline	Chlortetracycline
Concentration (parts per billion)			
0	0	0	–
100	18	42	–
150	53	89	–
200	80	99	–
250	94	100	49
300	98	100	99
600	–	–	–
Calculated 90/95	292	180	< 600 ppb

#### Cross-Reactivity

Concentration (parts per billion)	% Positive	
	Tetracycline	Doxycycline
0	0	0
30	0	0
60	0	0
150	100	100
300	100	100

The IDEXX SNAP Tetracycline Test Kit cross-reacts with the following drugs at the levels indicated: The IDEXX SNAP Tetracycline Test Kit does not cross-react with the following drugs at levels up to 100 ppb: Amoxicillin, Ampicillin, Ceftiofur, Cephapirin, Cloxacillin, Penicillin G, Sulfadiazine, Sulfanilamide, Sulfathiazole, Sulfamethazine, Sulfapyridine, Sulfadimethoxine, Gentamicin, Neomycin, Streptomycin, Ivermectin, Erythromycin, Pirimycin, Tilmicosin, Novobiocin, Fusosamide, Trichlormethazine, Thiabendazole, Chlorthiazide, Oxytocin, Phenylbutazone, Dexamethasone, Dipyrone, Flunixin, Enrofloxacin, PABA, Fluorfenicol, Minocycline.

### Selectivity

Sixty negative control milk samples were evaluated in an independent laboratory and none of these negative control samples tested positive with SNAP.

### Additional Test Information

Samples of this test kit model were independently evaluated by Antel Biosystems and were found to perform to the producer's specifications as stated in the test kit's descriptive insert.

### Sample Information

- Raw, commingled cow milk must be used.
- Samples must be refrigerated and testing must be completed within 72 hours of sampling.
- Thoroughly mix the sample before testing.
- For NCIMS purposes, milk samples cannot be frozen at any time. The testing of milk samples must be completed within 72 hours of sampling. Completed includes the time from initial test, any necessary retesting of the initial positive, confirmation and/or producer traceback testing.

### Checksets, Positive and Negative Controls in the United States for NCIMS Testing

Positive and negative controls must be run daily prior to testing samples, and with each new lot to verify the performance of the reagents and equipment. Read positive and negative checkset devices. Both devices must read within the limits on the checkset devices.

NOTE: Use the checkset with the black skirt for the SNAPshot Reader. Use the checkset with the blue skirt for the SNAPshot DSR Reader.

#### Negative Control (must be run daily):

Use tetracycline-negative raw milk; milk used for negative control purposes should read less than 0.95 on the IDEXX SNAPshot Reader. Negative control milk can be stored frozen for later use by aliquot and freezing at -15°C (5°F) or colder in a nonfrost-free freezer or in styrofoam container in a frost-free freezer. It must be used within 60 days. Thaw frozen milk at 0–4.5°C. Once thawed, mix thoroughly and do not use the aliquot if noticeable protein precipitation is present after thawing. Thawed negative control milk must be held at 0–4.5°C and can be used for 24 hours. Do not use beyond 24 hours. Milk controls may not be refrozen.

#### Positive Control (must be run daily):

Use the IDEXX SNAP ST Positive Control, part #98-27430-01, as directed in the descriptive insert that accompanies the product. The positive control should read greater than 1.2 on the IDEXX SNAPshot Reader. Do not dilute the positive control during confirmation testing. NOTE: If these conditions are not met, rerun the controls. If still not acceptable, contact IDEXX Technical Services.

### Technical assistance • [idexx.com/dairy](http://idexx.com/dairy)

US/Canada/Latinoamerica	+1 800 321 0207
EU	idexx.eu
中国	+010 84423133

日本	+021 6127 9528
Within AU	+1300 44 33 99
Within NZ	+0800 102 084

One IDEXX Drive • Westbrook, Maine 04092 USA

# IDEXX SNAP\* Tetracycline Test

The SNAP\* Tetracycline Test is validated to detect residues in raw, commingled cow milk. Contains biological components; in vitro use only.

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### Ensure sample and heater block are at correct temperature before use.

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


\*For NCIMS Appendix N Pilot Program please refer to the back page.


## 1 Transfer milk and swirl.

Transfiera la leche y agite.  
Transfira o leite e agite.  
Transférez le lait et mélangez.  
Przelej mleko i wymieszaj fiolkę.  
Milch transferieren und schwenken.  
Trasferire il latte e miscelare con un movimento circolare.  
倒入奶样并摇晃。  
生乳を移し、攪拌させます。


## 2 Incubate sample and SNAP test for 5 minutes

Incube la muestra y la prueba SNAP durante 5 minutos.  
Incube a amostra e o teste SNAP por 5 minutos.  
Laissez incuber l'échantillon et le test SNAP pendant 5 minutes.  
Poddawaj inkubacji próbkę i test SNAP przez 5 minut.  
Die Probe und den SNAP-Test 5 Minuten inkubieren.  
Incubare il campione e il test SNAP per 5 minuti.  
将样品和试剂盒放在加热器上孵育5分钟。  
サンプルと SNAP 検査品を 5 分間インキュベートします。


## 3 Pour, watch, then snap

Viértala, obsérvela y luego presione hasta que escuche el sonido SNAP.  
Despeje, observe e, então, faça o snap.  
Versez, observez, puis actionnez le SNAP.  
Przelej, obserwuj, a następnie docisnij.  
Abgießen, beobachten und dann herunterdrücken (Snap).  
Versare, osservare e premere finché si sente il suono di uno scatto ("snap" in inglese).  
倒入样品, 观察然后摁下。  
アクティベート サークル

### Snap when sample reaches activation circle.

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



OR  
0  
OU  
OU  
LUB  
ODER  
OPPURE  
或  
または

IDEXX SNAPshot\* DSR Reader

**Precautions:** Test milk within 3 days of collection; thoroughly mix milk before testing. Do not use abnormal-looking milk or expired tests; do not mix components from different kits. If the control spot fails to develop color, retest the sample.

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### It is recommended that positive samples be further analyzed.

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